


CANDO[®]

Develop strength and dexterity by progressing
from lighter to firmer colors



for a list of exercises
visit
www.FabEnt.com



Preferred colors[®] identify firmness

xx
soft

x
soft

soft

med

firm

x
firm



50-1010

Overdoor Shoulder Exerciser with Dual Pulleys

*Increase range-of-motion
Gain strength
Bracket fits over any door*

- helps increase and maintain range-of-motion and strength on all planes of shoulder movement
- stoppers prevent shoulder overextension
- adjustable length rope suitable for all exercise protocols and for any height user
- dual swivel pulleys assure smooth and easy operation
- bracket fits over any door

CANDO[®] 50-1010
exercise products




CORFLEX®

Designed for life.

Lace-Up Ankle with Stays

Recommended for preventing and preventing ankle sprains and instability

- > Designed to replicate ankle taping and fits easily inside most shoes
- > Features bilateral figure-8 strapping and removable plastic stays for added support
- > Manufactured from durable ballistic nylon for superior strength



Sizing Guide Measurement is circumference from base of heel around ankle							
2XS	XS	S	M	L	XL	2XL	3XL
9"-10"	10"-11"	11"-12"	12"-13"	13"-14"	14"-15"	15"-16"	16"+
(23-25cm)	(25-28cm)	(28-30cm)	(30-33cm)	(33-36cm)	(36-38cm)	(38-41cm)	(41+cm)



*Wear the products that
leading orthopedic and
sports medicine doctors
prescribe every day!*

800 426 7353 CORFLEX.COM

APPLICATION INSTRUCTIONS

1. Loosen laces and slide foot into brace.
2. Thread laces through eyelets, tighten and tie.
3. Grasp inside (medial) strap and wrap across top of foot, continue wrapping under heel and attach to side of the boot. Repeat for outside (lateral) strap.
4. Tighten appropriate strap to prevent inversion (outside strap) or eversion (inside strap) motion.
5. Wrap wide elastic strap over tied laces and secure.

Note: Brace should be snug but not hinder circulation. Stays can be removed to adjust support level.

CARE & CLEANING

Remove stays. Hand wash in warm water with mild detergent. Air dry.

EpiPoint®

ELLENBOGENSPANGE

ELBOW STRAP

BRACELET COUDIÈRE



**FIVE-POINT
PAD**

ENTLASTET DEN SEHNENANSATZ
RELIEVES THE TENDON INSERTION
SOULAGE LE POINT D'ATTACHE DU TENDON

FOUNDATION[®]

Shoulder Pulley



www.foundationrehabproducts.com